

## Density calibration files

The density calibration files used in The Image Processing Tool Kit and Fovea Pro are simple ascii text files containing the brightness to density calibration curve, and the values of the standards entered in the **IP\*Measure Global -> Calibrate Density** routine. The first line consists of the fixed labels "Brightness", "Density", "UserBrightness", and "CalibratedDensity", where the second label will be replaced by the name entered by the user. The first two columns contain brightness values from 0 to 255 and the corresponding density values. The third and fourth columns contain up to ten user-entered brightness and density standards. Columns are separated by a tab character.