

A COMPLETE INSTALLATION GUIDE TO

IHC PROFILER

CONTENTS

- I. ImageJ
- II. IHC Profiler
- III. Tips and Troubleshooting

I. ImageJ

- Download ImageJ from <http://rsb.info.nih.gov/ij/download.html>
- After downloading, open the installation file and follow the onscreen instructions and install the software respectively.


II. IHC Profiler – An Installation Guide

- Extract/Unzip the **IHC_Profiler.zip** file containing the **IHC Profiler** folder and a **IHC_Profiler.txt** macro file.
- Copy *both* the files (IHC Profiler folder and IHC_Profiler text file) and paste it in the following location or similar depending on where ImageJ has been installed:
 - ❖ My Computer > C: > Program Files > ImageJ > plugins
- Additionally, copy *just* the **IHC_Profiler.txt** file and paste it into the following location or similar depending on where ImageJ has been installed:
 - ❖ My Computer > C: > Program Files > ImageJ > macros
- Restart ImageJ (if already open).


III. Tips and Troubleshooting

- We recommend using ImageJ version 1.45s or higher.
- Depending upon the system specification, occurrence of false profile plots might be an issue due to improper or low buffering rate of the computer system. We suggest refreshing the screen at regular intervals and opening one image at a time for analysis for better and accurate score profile results.
- If you continue to get results which you think may be false based on your visual perception, kindly reset clipboard and buffer as follows:

Plugins > Utility > Reset



Click on **Reset...**



Using the drop down menu

Select Reset: **Clipboard** and Click **OK**

Select Reset: **Undo Buffer** and Click **OK**

- If the problem persists, you may try restarting your computer.
- For better performance results update your system to Java 6 or higher.